



## Published Applications

### IMMUNOLOGY

· Activates T cells from human, mouse, or rat sources, in combination with PMA, to express cytokines including IL-17, IL-4, IL-10, and IL-2 (Caraher et al.; Harrington et al.; Parrish-Novak et al.).

## References

- Caraher EM et al. (2000) Flow cytometric analysis of intracellular IFN-gamma, IL-4 and IL-10 in CD3(+)4(+) T-cells from rat spleen. *J Immunol Methods* 244(1-2): 29–40.
- Elzi DJ et al. (2001) Ionomycin causes activation of p38 and p42/44 mitogen-activated protein kinases in human neutrophils. *Am J Physiol Cell Physiol* 281(1): C350–60.
- Harrington LE et al. (2005) Interleukin 17-producing CD4+ effector T cells develop via a lineage distinct from the T helper type 1 and 2 lineages. *Nature Immunol* 6(11): 1123–32.
- Liu C & Hermann TE. (1978) Characterization of ionomycin as a calcium ionophore. *J Biol Chem* 253(17): 5892–4.
- Liu WC et al. (1978) Ionomycin, a new polyether antibiotic. *J Antibiot (Tokyo)* 31(9): 815–9.
- Morgan AJ & Jacob R. (1994) Ionomycin enhances Ca<sup>2+</sup> influx by stimulating store-regulated cation entry and not by a direct action at the plasma membrane. *Biochem J* 300(Pt 3): 665–72.
- Parrish-Novak J et al. (2000) Interleukin 21 and its receptor are involved in NK cell expansion and regulation of lymphocyte function. *Nature* 408: 57–63.
- Yoshida S & Plant S. (1992) Mechanism of release of Ca<sup>2+</sup> from intracellular stores in response to ionomycin in oocytes of the frog *Xenopus laevis*. *J Physiol* 458: 307–18.
- Zhang LH et al. (1997) Antiproliferative and immunosuppressive properties of microcolin A, a marine-derived lipopeptide. *Life Sci* 60(10): 751–62.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

**This product is hazardous. Please refer to the Safety Data Sheet (SDS).**

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.