

## Small Molecules

D4476

Activin/NODAL/TGF- $\beta$  pathway inhibitor; Inhibits CK1, ALK5 and p38 MAPK

Catalog # 74012  
74014

5 mg  
10 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

D4476 is a selective, potent, and cell-permeant inhibitor of casein kinase (CK) 1 (Rena et al.). It inhibits CK1 and CK1 $\delta$  from *S. pombe* (IC<sub>50</sub> = 200 and 300 nM, respectively; Bain et al.; Rena et al.). It also inhibits transforming growth factor  $\beta$ 1 (TGF- $\beta$ 1) type I receptor (ALK5) with IC<sub>50</sub> = 500 nM (Rena et al.).

Molecular Name:	D4476
Alternative Names:	Casein Kinase I Inhibitor
CAS Number:	301836-43-1
Chemical Formula:	C <sub>23</sub> H <sub>18</sub> N <sub>4</sub> O <sub>3</sub>
Molecular Weight:	398.4 g/mol
Purity:	≥ 98%
Chemical Name:	4-[4-(2,3-dihydro-1,4-benzodioxin-6-yl)-5-(2-pyridinyl)1H-imidazol-2-yl]-benzamide
Structure:	



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 75 mM · Absolute ethanol ≤ 25 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 251 $\mu$ L of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

## Published Applications

### REPROGRAMMING

- Along with Forskolin (Catalog #72112) and 2-methyl-5-hydroxytryptamine, reprograms mouse somatic cells into pluripotent cells by acting as a chemical substitute for OCT4 (Hou et al.).
- Reprograms mouse somatic cells into pluripotent cells by acting as a chemical substitute for OCT4 (Hou et al.).

### CANCER RESEARCH

- Induces apoptosis of multiple myeloma cell lines (Hu et al.).
- Exhibits anti-leukemic activity by selectively killing leukemia stem cells over normal hematopoietic stem and progenitor cells (Järås et al.).

## References

- Bain J et al. (2003) The specificities of protein kinase inhibitors: an update. *Biochem J* 371(1): 199–204.
- Hou P et al. (2013) Pluripotent stem cells induced from mouse somatic cells by small-molecule compounds. *Science* 341(6146): 651–4.
- Hu Y et al. (2015) CSNK1 $\alpha$ 1 mediates malignant plasma cell survival. *Leukemia* 29(2): 474–82.
- Järås M et al. (2014) Csnk1a1 inhibition has p53-dependent therapeutic efficacy in acute myeloid leukemia. *J Exp Med* 211(4): 605–12.
- Rena G et al. (2004) D4476, a cell-permeant inhibitor of CK1, suppresses the site-specific phosphorylation and nuclear exclusion of FOXO1a. *EMBO Rep* 5(1): 60–5.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.