

## Small Molecules

### Geldanamycin

Inhibits Hsp90

Catalog # 74032  
74034

1 mg  
5 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM


FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Geldanamycin is a benzoquinone ansamycin antibiotic that selectively inhibits heat shock protein 90 (Hsp90) (Schulte et al.; Whitesell et al.). It inhibits the ATPase activity of Hsp90 by binding with high affinity to the N-terminal ATP-binding site ( $K_d = 1.2 \mu\text{M}$ ; Roe et al.). Geldanamycin also binds to the glucose-regulated protein GRP94, an endoplasmic reticulum protein with homology to Hsp90. This interaction destabilizes its complex with p185-ErbB2, leading to degradation of ErbB2 protein, which is overexpressed in aggressive breast cancers (Castagnola et al.; Chavany et al.).

Molecular Name:	Geldanamycin
Alternative Names:	NSC 122750
CAS Number:	30562-34-6
Chemical Formula:	$\text{C}_{29}\text{H}_{40}\text{N}_2\text{O}_9$
Molecular Weight:	560.6 g/mol
Purity:	$\geq 98\%$
Chemical Name:	(6-hydroxy-5,11,21-trimethoxy-3,7,9,15-tetramethyl-16,20,22-trioxo-17-azabicyclo[16.3.1]docosa-1(21),8,12,14,18-pentaen-10-yl) carbamate

Structure:



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at $-20^\circ\text{C}$ as supplied. Protect product from prolonged exposure to light. For long-term storage store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO $\leq 15 \text{ mM}$ For example, to prepare a 1 mM stock solution in DMSO, resuspend 1 mg in 1.8 mL of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at  $-20^\circ\text{C}$ . Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

## Published Applications

### CANCER RESEARCH

- Inhibits growth of 60 human tumor cell lines (Supko et al.).

### DISEASE MODELING

- Regulates innate immune response as demonstrated by inhibition of inflammasome activity in a mouse gout model (Mayor et al.).
- Inhibits replication of herpes simplex virus (HSV) 1 and HSV-2 (Li et al.).

## References

- Castagnola P et al. (2016) Identification of an HSP90 modulated multi-step process for ERBB2 degradation in breast cancer cells. *Oncotarget* 7(51): 85411–29.
- Chavany C et al. (1996) p185erbB2 binds to GRP94 in vivo. Dissociation of the p185erbB2/GRP94 heterocomplex by benzoquinone ansamycins precedes depletion of p185erbB2. *J Biol Chem* 271(9): 4974–7.
- Li Y-H et al. (2012) Geldanamycin, a ligand of heat shock protein 90, inhibits herpes simplex virus type 2 replication both in vitro and in vivo. *J Antibiot* 65(10): 509–12.
- Mayor A et al. (2007) A crucial function of SGT1 and HSP90 in inflammasome activity links mammalian and plant innate immune responses. *Nat Immunol* 8(5): 497–503.
- Roe SM et al. (1999) Structural basis for inhibition of the Hsp90 molecular chaperone by the antitumor antibiotics radicicol and geldanamycin. *J Med Chem* 42(2): 260–6.
- Schulte TW et al. (1995) Disruption of the Raf-1-Hsp90 molecular complex results in destabilization of Raf-1 and loss of Raf-1-Ras association. *J Biol Chem* 270(41): 24585–8.
- Supko JG et al. (1995) Preclinical pharmacologic evaluation of geldanamycin as an antitumor agent. *Cancer Chemother Pharmacol* 36(4): 305–15.
- Whitesell L et al. (1994) Inhibition of heat shock protein HSP90-pp60v-src heteroprotein complex formation by benzoquinone ansamycins: essential role for stress proteins in oncogenic transformation. *Proc Natl Acad Sci USA* 91(18): 8324–8.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.