

# Round-Bottom Polystyrene Tubes without Caps, 5 mL

Polystyrene tubes without caps, 1000/case

Catalog # 38055

1000 Tubes



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Falcon® Round-Bottom Polystyrene Tubes, 5 mL are sterile, high-quality tubes ideal for a wide variety of sample preparation, sample storage, and sample centrifugation tasks. They are required for EasySep™ cell isolation when using EasySep™ Magnet (Catalog #18000) or EasyEights™ EasySep™ Magnet (5 mL side; Catalog #18103), and can be used for flow cytometry and fluorescence-activated cell sorting (FACS). Tubes are RNase-/DNase-free.

500 Tubes with caps; 25 Tubes/pack (Corning Catalog #352058)

1000 Tubes with caps; 125 Tubes/pack (Corning Catalog #352054)

1000 Tubes without caps; 125 Tubes/pack (Corning Catalog #352052)

### APPLICATIONS

- Cell separation
- Cell culture
- Sample storage
- Liquid handling

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Clear and transparent polystyrene
- Biologically inert, RNase-/DNase-free, and non-pyrogenic

### DIMENSIONS

- 12 x 75 mm (outer diameter x height)

### NOMINAL VOLUME

- 5 mL
- No graduations

### CENTRIFUGAL FORCE

- Up to 1400 x g

### STERILITY

- Gamma irradiation

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasyEights, and EasySep are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.