

# Cell Scraper, Small

Small cell scraper, 24 cm

Catalog # 38065

100 Scrapers



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Cell Scraper is sterile, individually wrapped, non-pyrogenic, disposable, and have a polystyrene handle and a low-density polyethylene blade. They are useful for harvesting cells or cell lysates and are designed for use with Tissue Culture-Treated Dishes (Catalog #38046) and Cell Culture Flasks (Catalog #38014).

Small - 100, individually wrapped (Trueline Catalog #TR9000 and equivalent to Corning Catalog #3010/353086)

Large - 100, individually wrapped (Trueline Catalog #TR9001 and equivalent to Corning Catalog #3011/353087)

### APPLICATIONS

- Cell Culture

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** For product expiry date, contact [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Features

### MATERIAL

- Handle: green polystyrene
- Blade: white low-density polyethylene
- Non-pyrogenic

### DIMENSIONS

- 1.8 cm blade, 24 cm handle

### STERILITY

- Irradiation
- Individually wrapped

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.