

Cryogenic Vial Cap Inserts, Assorted Colors



Colored cryogenic vial cap inserts in white, blue, red, green, and yellow for sample color-coding

Catalog # 38081

500 Inserts

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Corning® Cryogenic Vial Cap Inserts come in white, blue, red, green, and yellow in a resealable bag. Cap inserts are useful for color-coding vials for easy sample identification. Non-sterile polypropylene inserts are designed to fit most brands of cryogenic vials (e.g. Catalog #38047, 38048, 38049, or 38053).

500; 100 of each color (Corning Catalog #430499)

APPLICATIONS

- Cryopreservation
- Cell storage

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For product expiry date, contact techsupport@stemcell.com.

Features

- FIT
- Fits all cryogenic vials carried by STEMCELL and most standard cryogenic vials

DIMENSIONS

- Diameter 9.6 mm
- Length 5.1 mm

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Corning is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.