

Cell Spreader, T Form

T-shaped cell spreader

Catalog # 38082

3500 Each



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Cell Spreader (also known as Gosselin™ T Form green polystyrene cell spreader) is sterile and has rounded ends to enable smooth contact with growth surfaces. The T-shape ensures even pressure and is designed for spreading and dispersing liquids across surfaces such as agar. Packaged in resealable bags.

3500 Spreaders; 25/Bag (Corning Catalog #ETAT-04)

APPLICATIONS

- Cell culture

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Green polystyrene

DIMENSIONS

- Length: 150 mm
- Shape: T-form

STERILITY

- Sterile

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Gosselin is a trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.