

Polystyrene Storage Bottle, 125 mL



Polystyrene bottle, 125 mL

Scientists Helping Scientists™ | WWW.STEMCELL.COM

Catalog # 38083

24 Bottles

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Costar® Polystyrene Storage Bottle comes sterile with a high-density polyethylene plug seal screw cap that provides an airtight seal to reduce risk of contamination. Smooth-sided bottles with a wide mouth (45 mm diameter) are made of high-clarity, non-pyrogenic polystyrene and are ideal for storing media, buffers, and other aqueous solutions.

125 mL - 24 Bottles, individually wrapped (Corning Catalog #8388)

250 mL - 12 Bottles, individually wrapped (Corning Catalog #8390)

APPLICATIONS

- Liquid storage

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Bottle: polystyrene
- Bottle shape: cylindrical
- Cap: red high-density polyethylene
- Non-pyrogenic

DIMENSIONS

- Height: 103.4 mm (with cap)
- Neck: 38.6 x 44.2 mm (inner diameter x outer diameter with threads)
- Cap: 45 mm

NOMINAL VOLUME

- 125 mL
- Graduated markings

STERILITY

- Gamma irradiation
- Individually wrapped

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Costar is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.