

Polyethylene Storage Bottle, 125 mL



Polyethylene bottle, 125 mL

Catalog # 38085

336 Bottles

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Polyethylene Storage Bottle (also known as Gosselin™ 125 mL storage bottle) is made of polyethylene terephthalate (PET), has a transparency equivalent to glass, and is sterile, lightweight, and break-resistant. These bottles are ideal for storing media, buffers, and other aqueous solutions. Each bottle comes with a leak-proof, tamper-evident, high-density polyethylene screw cap.

336 Bottles (Corning Catalog #P125B-50)

APPLICATIONS

- Liquid storage

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Bottle: optically clear polyethylene terephthalate (PET)
- Bottle shape: octagonal
- Cap: natural color high-density polyethylene

DIMENSIONS

- Height : 105.0 mm (with cap)
- Cross section: 53 mm
- Cap: 31.7 mm

NOMINAL VOLUME

- 125 mL
- Graduations every 25 mL

STERILITY

- Sterile

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Gosselin is a trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.