

# Transfer Pipettes, 1 mL



## Transfer pipettes, 1 mL

Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

Catalog # 38086

2000 Pipettes

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Transfer Pipettes by Gosselin™ are sterile, individually wrapped, disposable pipettes designed for quick and dropwise liquid transfer. Built-in pipette bulb and metric graduations on each transfer pipette enables efficient handling of liquids and suspensions without the need for any other tools. These polyethylene pipettes have reduced cell and protein binding compared to glassware, and are sterilized with an ethylene oxide treatment.

1 mL - 2000 Pipettes, individually wrapped (Corning Catalog #PAS1A-01)

3 mL - 2000 Pipettes, individually wrapped (Corning Catalog #PAS3A-01)

### APPLICATIONS

- Liquid handling

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Natural color polyethylene with built-in bulb

### DIMENSIONS

- Length: 150 mm
- Metric graduations

### CAPACITY

- 1 mL

### STERILITY

- Ethylene oxide treatment
- Individually wrapped

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Gosselin is a trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.