

# Falcon® Serological Pipettes, 1 mL

**Sterile polystyrene pipettes for cell culture**

Catalog # 38001

1000 Pipettes



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Falcon® Serological Pipettes, 1 mL can be used for tissue culture or with any non-toxic solution that requires aliquoting or aspirating/dispensing in a sterile environment. These pipettes have been tested for dispensing and aspirating mTeSR™1 (Catalog #85850) and TeSR™2 (Catalog #05860) media. They have also been used in cell separation with the EasySep™ Magnet (Catalog #18000) and EasyEight™ EasySep™ Magnet (Catalog #18103).

1000 Pipettes; 100 Pipettes/box, individually wrapped (Corning Catalog #357521)

### APPLICATIONS

- Cell separation
- Tissue culture
- Aliquoting or aspirating/dispensing in a sterile environment

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Pipette: clear and transparent polystyrene
- Plug: polyester fiber
- Color-coded yellow
- Non-pyrogenic, low-endotoxin

### CAPACITY

- 1 mL
- Increments of 0.01 mL with reverse graduations
- Overfill is 0.04 mL
- Reliable and accurate to  $\pm 3\%$  at full volume

### STERILITY

- Gamma irradiation

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasyEights, and EasySep are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. mTeSR and TeSR are trademarks of WARF. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.