

# Costar® 6-Well Flat-Bottom Plate, Tissue Culture-Treated

Clear polystyrene flat-bottom, tissue culture-treated multiwell plate with lid

Catalog # 38015      50 Plates



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

The Costar® 6-Well Flat-Bottom Plate, Tissue Culture-Treated is recommended for use with NeuroCult™ NS-A Proliferation Kit (Human) (Catalog #05751). This plate can also be used in cell culture or in any application where flat-bottom plates with a hydrophilic and negatively charged surface are needed.

50 Plates, individually wrapped (Corning Catalog #3516)

### APPLICATIONS

- Cell culture

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Plate and lid: clear and transparent polystyrene
- Hydrophilic and negatively charged
- Non-pyrogenic

### DIMENSIONS

- Conforms to standard plate footprint and dimensions

### CAPACITY

- 16.8 mL
- Recommended working volume of 1.9 - 2.9 mL per well

### STERILITY

- Gamma irradiation

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For information on chemical resistance or thermoplastic property, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and NeuroCult are trademarks of STEMCELL Technologies Canada Inc. Costar is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.