

Falcon® 6-Well Flat-Bottom Plate, Tissue Culture-Treated

Clear polystyrene, flat-bottom, tissue culture-treated multiwell cell culture plate with lid

Catalog # 38016

50 Plates



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Falcon® 6-Well Flat-Bottom Plate, Tissue Culture-Treated is recommended for use with mTeSR™1 (Catalog #85850). This plate can also be used to culture anchorage-dependent mammalian cells, including stromal cells, and in any application where aqueous reagents are needed for tissue culture.

50 Plates, individually wrapped (Corning Catalog #353046)

APPLICATIONS

- Cell culture

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For product expiry date, contact techsupport@stemcell.com.

Features

MATERIALS

- Plate and lid: clear and transparent polystyrene
- Hydrophilic and negatively charged
- Non-pyrogenic

DIMENSIONS

- Conforms to standard plate footprint and dimensions

CAPACITY

- 15.5 mL
- Recommended working volume of 1.9 - 2.9 mL per well

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. mTeSR is a trademark of WARF. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.