

Falcon® 96-Well Flat-Bottom Microplate, Tissue Culture-Treated

Clear polystyrene flat-bottom, tissue culture-treated multiwell plate
with lid

Catalog # 38022

50 Plates



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

The Falcon® 96-Well Flat-Bottom Microplate, Tissue Culture-Treated has been tested for use with StemSpan™ SFEM II (Catalog #09605), StemSpan™ SFEM (Catalog #09600), and StemSpan™ H3000 (Catalog #09800). This plate is also recommended for monolayer cell culture and immunoassays, or in any application where clear, flat-bottom, polystyrene plates are used, such as cell imaging, cell culture with binding medium (> 10 kDa), or cell culture with large biomolecules that possess ionic groups and/or hydrophobic regions.

50 Plates, individually wrapped (Corning Catalog #353072)

APPLICATIONS

- Cell culture
- Immunoassays
- Cell imaging

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For product expiry date, please contact techsupport@stemcell.com.

Features

MATERIALS

- Plate and lid: clear and transparent polystyrene
- Hydrophilic and negatively charged
- Non-pyrogenic

DIMENSIONS

- Conforms to standard plate footprint and dimensions

CAPACITY

- 0.37 mL per well
- Recommended working volume of 0.04 - 0.27 mL

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For information on chemical resistance or thermoplastic property, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and StemSpan are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.