

# ELISA Diluent



**Diluent for ELISA samples, standards, and detection antibodies**

Catalog # 01916

2 x 120 mL

Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

ELISA Diluent is intended for dilution of samples, standards, and detection antibodies in ELISA assays. By preventing the cross-linking of heterophilic antibodies, ELISA Diluent minimizes interference and false-positive signals in ELISA reactions. ELISA Diluent is required for use with the Human IgE ELISA Antibody Pair Kit (Catalog #01993) when serum or plasma samples are used.

NOTE: Do not use ELISA Diluent for dilution of streptavidin conjugates.

## Properties

- Storage:** Store at 2 - 8°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- Protein-containing buffer
  - 0.002% Kathon CG

## Handling / Directions For Use

Refer to the Product Information Sheet for the applicable ELISA kit.

### RELATED PRODUCTS

For a complete list of ELISA kits and related products from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.