

Corning® Filtered Pipette Tips, 1000 µL



Sterile polypropylene filtered pipette tips

Scientists Helping Scientists™ | WWW.STEMCELL.COM

Catalog # 38031

1000 Tips

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Corning® Filtered Pipette Tips are ideal for precise micropipetting tasks in a sterile environment. These pipette tips fit most major-brand micropipettes and are filtered with a non-reactive, hydrophobic material to prevent pipette cross-contamination. Tips are pre-racked in 10 polypropylene containers per case.

2 µL - 960 Tips in 10 Racks (Corning Catalog #4801)

10 µL - 960 Tips in 10 Racks (Corning Catalog #4807 and equivalent to Corning Catalog #GEN3-10-L-R-S/TF-400)

30 µL - 960 Tips in 10 Racks (Corning Catalog #4821)

200 µL - 960 Tips in 10 Racks (Corning Catalog #4810 and equivalent to Corning Catalog #GEN3-200-L-R-S)

1000 µL - 1000 Tips in 10 Racks (Corning Catalog #4809 and equivalent to Corning Catalog #GEN3-1000-L-R-S/TF-1000-R-S)

APPLICATIONS

- Liquid handling
- RNA/DNA purification
- PCR
- Molecular biology

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Tip: polypropylene
- Filter: hydrophobic, non-reactive fiber
- Pre-racked in orange polypropylene rack
- RNase/DNase-free and non-pyrogenic

DIMENSIONS

- Total tip length of 7.234 cm
- Internal diameter of 0.750 cm at large opening
- Internal diameter of 0.091 cm at small opening

CAPACITY

- 100 - 1000 µL
- Graduations at 250 µL, 500 µL, and 1000 µL

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Corning is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.