

48-Well Flat-Bottom Plate, Non-Treated



Sterile, flat-bottom, non-treated 48-well plate with lid

Scientists Helping Scientists™ | WWW.STEMCELL.COM

Catalog # 38043

50 Plates

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Flat-Bottom Plates, Non-Treated are ideal for a variety of cell culture applications. The plates have a clear, untreated, hydrophobic polystyrene surface and are sterile and ready for use in an aseptic environment. The included lid is non-reversible and is designed to prevent cross-condensation among wells. These plates are stackable for easy handling and storage.

6-Well - 100 Plates; 5 Plates/bag (Corning Catalog #3736 and equivalent to Corning Catalog #351146)

12-Well - 100 Plates; 5 Plates/bag (Corning Catalog #3737 and equivalent to Corning Catalog #351143)

24-Well - 100 Plates; 5 Plates/bag (Corning Catalog #3738 and equivalent to Corning Catalog #351147)

48-Well - 50 Plates, individually wrapped (Corning Catalog #351178)

96-Well - 50 Plates, individually wrapped (Corning Catalog #351172)

APPLICATIONS

- Cell culture

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For product expiry date, contact techsupport@stemcell.com.

Features

MATERIALS

- Plate and lid: clear and transparent polystyrene
- RNase/DNase-free and non-pyrogenic

DIMENSIONS

- 127.6 x 85.5 mm (length x width)
- Conforms to standard plate footprint and dimensions
- Stackable for easy handling and storage

CAPACITY

- 1.4 mL per well
- Recommended working volume of 0.19 - 0.285 mL per well
- Approximate cell growth area of 0.75 cm² per well

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.