

# 60 mm Dish, Tissue Culture-Treated



**Sterile 60 mm tissue culture-treated dish**

Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

Catalog # 38051

500 Dishes

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Tissue Culture-Treated Dishes are ideal for attachment-dependent cell and tissue culture. These polystyrene dishes are transparent for superior optical clarity and are stackable for easy handling and storage. The dishes have been sterilized and are vented to allow for gas exchange.

35 mm - 500 Dishes; 20 Dishes/bag (Trueline Catalog #TR4000 and equivalent to Corning Catalog #353001/430165)

60 mm - 500 Dishes; 20 Dishes/bag (Trueline Catalog #TR4001 and equivalent to Corning Catalog #353002/430166)

100 mm - 500 Dishes; 20 Dishes/bag (Corning Catalog #430167 and equivalent to Corning Catalog #430293/353003)

### APPLICATIONS

- Cell culture
- Microscopy

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** For product expiry date, contact [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Features

### MATERIALS

- Dish and lid: clear and transparent polystyrene
- Tissue culture-treated for optimal cell attachment
- Non-pyrogenic

### DIMENSIONS

- 60 x 15 mm (diameter x height)
- Stackable for easy handling and storage

### CAPACITY

- Recommended working volume of 4.2 - 6.3 mL
- Approximate cell growth area of 21 cm<sup>2</sup>

### STERILITY

- Gamma irradiation

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.