

# Corning® Cryogenic Vials with Orange Caps, 2 mL

**Sterile, self-standing, round-bottom, internally-threaded, polypropylene cryogenic vial with orange cap**

Catalog # 38053

500 Vials



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Corning® Cryogenic Vials with Orange Caps, 2 mL are ideal for reliable sample cryostorage.

500 Vials; 50 Vials/pack (Corning Catalog #430488)

### APPLICATIONS

- Cryopreservation
- Cell storage

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIALS

- Tube: clear and transparent polypropylene, round-bottom, self-standing
- Cap: orange polypropylene, internally threaded
- Washer: silicone
- Withstands temperatures as low as -196°C (in gas/vapor phase)
- RNase/DNase-free and non-pyrogenic

### DIMENSIONS

- Length - Vial with cap: 48.5 mm
- Diameter - Vial: 12.7 mm; Cap: 11.9 mm

### NOMINAL VOLUME

- 2 mL
- Graduations at 1.0 mL and 1.8 mL

### STERILITY

- Gamma irradiation

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Corning is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.