

Cytokines

Human Recombinant VEGF-D



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Vascular endothelial growth factor D

Catalog # 78203
78203.1

10 µg
50 µg

Product Description

Vascular endothelial growth factor D (VEGF-D) is a member of the VEGF/platelet-derived growth factor (PDGF) family of proteins. VEGF-D is a potent angiogenic factor and promotes lymphangiogenesis, endothelial cell growth and survival, and can affect blood vessel permeability. VEGF-D is expressed in the lung, heart, small intestine, fetal lung, and at lower levels in the pancreas, colon, and skeletal muscle (Otrock et al.; Roy et al.; Stacker et al.; Yamada et al.). VEGF-D is a ligand for VEGF receptors 2 (VEGFR-2 [Flk1]) and 3 (VEGFR-3 [Flt4]) (Baldwin et al.). VEGFR-3 is highly expressed in lymphatic endothelial cells and is essential for their growth and differentiation (Otrock et al.; Roy et al.). Binding of VEGF-D to neuropilins contributes to VEGFR-3 signaling during lymphangiogenesis, whereas binding to integrin $\alpha 9 \beta 1$ promotes endothelial cell adhesion and migration (Roy et al.; Otrock et al.). During embryogenesis, VEGF-D also plays a role in the formation of the venous and lymphatic systems.

Product Information

Alternative Names: c-Fos induced growth factor, FIGF, Vascular endothelial growth factor D
Accession Number: O43915
Amino Acid Sequence: FAATFYDIET LKVIDEEWQR TQCSPRETCV EVASELGKST NTFFKPPCVN VFRCGGCCNE ESLICMNTST SYISKQLFEI SVPLTSVPEL VPKVANHTG CKCLPTAPRH PYSIIRR
Predicted Molecular Mass: 13.1 kDa
Species: Human
Cross Reactivity: Reported to be species-specific
Formulation: Lyophilized after dialysis against phosphate-buffered saline.
Source: CHO


Specifications

Activity: The specific activity is $\geq 1.0 \times 10^3$ units/mg ($EC_{50} \leq 1000$ ng/mL) as determined by a cell proliferation assay using human umbilical vein endothelial cells (HUVECs).
Purity: $\geq 95\%$
Endotoxin Level: Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is ≤ 0.2 EU/µg protein.

Preparation and Storage

Storage: Store at -80°C .
Stability: Stable as supplied for 12 months from date of receipt.
Preparation: Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex. As a general guide, do not store at $2 - 8^\circ\text{C}$ for more than 1 week or at -20°C for more than 2 months. Avoid repeated freeze-thaw cycles.

Data



(A) The biological activity of Human Recombinant VEGF-D was tested by its ability to promote the proliferation of HUVECs. Cell proliferation was measured using a fluorometric assay method. The EC₅₀ is defined as the effective concentration of the growth factor at which cell proliferation is at 50% of maximum. The EC₅₀ in the example above is less than 1000 ng/mL.

(B) 2 µg of Human Recombinant VEGF-D was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Human Recombinant VEGF-D has a predicted molecular mass of 13.1 kDa.

Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/cytokines or contact us at techsupport@stemcell.com.

References

- Baldwin ME et al. (2001) The specificity of receptor binding by vascular endothelial growth factor-d is different in mouse and man. *J Biol Chem* 276(22): 19166–71.
- Otrock ZK et al. (2007) Vascular endothelial growth factor family of ligands and receptors: review. *Blood Cells Mol Dis* 38(3): 258–68.
- Roy H et al. (2006) Biology of vascular endothelial growth factors. *FEBS Lett* 580(12): 2879–87.
- Stacker SA et al. (1999) Biosynthesis of vascular endothelial growth factor-D involves proteolytic processing which generates non-covalent homodimers. *J Biol Chem* 274(45): 32127–36.
- Yamada Y et al. (1997) Molecular cloning of a novel vascular endothelial growth factor, VEGF-D. *Genomics* 42(3): 483–8.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.