

# NeuroCult™ SM1 Without Antioxidants



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog # 05732

10 mL

## Product Description

NeuroCult™ SM1 Without Antioxidants is an optimized serum-free culture supplement, based on the published B27 formulation (Brewer et al.), with antioxidants (tocopherol, tocopherol acetate, superoxide dismutase, catalase, and glutathione) removed. NeuroCult™ SM1 Without Antioxidants can be used in protocols where antioxidants may interfere with oxidative stress, or similar assays used in neurodegenerative research. In combination with a basal medium of choice, NeuroCult™ SM1 Without Antioxidants can be used for neural progenitor cell expansion and the study of neural stem and progenitor cell differentiation.

- Versatile cell culture supplement
- Optimized, serum-free formulation
- Raw materials rigorously screened

## Properties

<b>Storage:</b>	Store at -20°C.
<b>Shelf Life:</b>	Stable until expiry date (EXP) on label.
<b>Contains:</b>	<ul style="list-style-type: none"><li>• Insulin</li><li>• Other ingredients</li></ul>

This product contains material derived from human plasma. Donors have been tested and found negative for HIV-1 and -2, hepatitis B, and hepatitis C prior to donation. However, this product should be considered potentially infectious and treated in accordance with universal handling precautions.

Lot-to-lot variability in color of this product may be expected. This will not affect performance.

## Handling / Directions For Use

1. Thaw bottle of NeuroCult™ SM1 Without Antioxidants at room temperature (15 - 25°C) for 1 hour or at 2 - 8°C overnight.  
NOTE: If not used immediately, aliquot and store at -20°C. Do not exceed the expiry date (EXP) as indicated on label.
2. Add supplement to basal medium (e.g. DMEM/F-12 with 15 mM HEPES; Catalog #36254) at a 1 in 50 dilution. Mix thoroughly.  
NOTE: If not used immediately, store at 2 - 8°C for up to 1 month.
3. Use as directed in the protocol of choice.

## Notes and Tips

For related products, including specialized culture media, cytokines, dissociation reagents, and cultureware, visit [www.stemcell.com/NSPCworkflow](http://www.stemcell.com/NSPCworkflow) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

- Brewer GJ et al. (1993) Optimized survival of hippocampal neurons in B27-supplemented Neurobasal, a new serum-free medium combination. *J Neurosci Res* 35(5): 567–76.
- Brewer GJ & Cotman CW. (1989) Survival and growth of hippocampal neurons in defined medium at low density: advantages of a sandwich culture technique or low oxygen. *Brain Res* 494(1): 65–74.
- Chen Y et al. (2008) NS21: re-defined and modified supplement B27 for neuronal cultures. *J Neurosci Methods* 171(2): 239–47.
- O'Connor TJ et al. (1998) Isolation and propagation of stem cells from various regions of the embryonic mammalian central nervous system. In: Celis JE (Ed.) *Cell Biology. A Laboratory Handbook* (p. 149). London, UK: Academic Press.
- Perry SW et al. (2004) Antioxidants are required during the early critical period, but not later, for neuronal survival. *J Neurosci Res* 78(4): 485–92.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and NeuroCult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.