

# Axygen® PCR Domed Strip Caps, Clear, 8 Caps/Strip

Non-sterile, clear, polypropylene domed PCR strip caps

Catalog # 38095

1250 Strips



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Axygen® PCR Strip Caps are RNase/DNase-free, non-pyrogenic, non-sterile, and are ideal for polymerase chain reaction (PCR) and real-time PCR (qPCR). The cap design allows a tight seal, reduces sample evaporation, and is compatible with Axxygen® PCR Strip Tubes (Catalog #38094) and Axxygen® PCR Microplates (Catalog #38098). Axxygen® PCR Strip Caps are available in domed, flat, or flat/ultra-clear designs. The flat cap allows easy labeling of samples, the domed cap increases efficiency of heat transfer, and the flat/ultra-clear cap allows for 75% higher light transmittance.

Catalog #38095, Domed - 1250 Strips; 125 Strips/pack (Corning Catalog #PCR-02CP-C)

Catalog #38096, Flat - 1250 Strips; 125 Strips/pack (Corning Catalog #PCR-02-FCP-C)

Catalog #38097, Flat, Ultra-Clear - 1250 Strips; 125 Strips/pack (Corning Catalog #PCR-2CP-RT-C)

### APPLICATIONS

- PCR
- Real-time PCR (qPCR)

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on manufacturer label.

## Features

### MATERIAL

- Clear polypropylene
- RNase/DNase-free and non-pyrogenic
- Domed cap

### STERILITY

- Non-sterile

## Handling / Directions For Use

Use as directed in protocol of choice.

## Notes and Tips

For chemical resistance or thermoplastic property information, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Axxygen is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.