

Small Molecules

(Z)-4-Hydroxytamoxifen

Estrogen signaling pathway modulator; Modulates selective estrogen receptor and estrogen-related receptor

Catalog #	74052	5 mg
	74054	10 mg
	74056	25 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

(Z)-4-Hydroxytamoxifen is a cell-permeable selective estrogen receptor modulator (SERM) that binds to estrogen receptors (ER) and estrogen-related receptors (ERR) with both estrogenic and anti-estrogenic effects (Poulin et al.).

Molecular Name:	(Z)-4-Hydroxytamoxifen
Alternative Names:	trans-4-Hydroxytamoxifen; 4-OH-TAM; ICI 79280
CAS Number:	68047-06-3
Chemical Formula:	C ₂₆ H ₂₉ NO ₂
Molecular Weight:	387.5 g/mol
Purity:	≥ 98%
Chemical Name:	4-[(1Z)-1-[4-[2-(dimethylamino)ethoxy]phenyl]-2-phenylbut-1-en-1-yl]phenol
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 5.1 mM · Absolute ethanol ≤ 50 mM · DMF ≤ 50 mM For example, to prepare a 10 mM stock solution in absolute ethanol, resuspend 1 mg in 258 µL of absolute ethanol.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in absolute ethanol at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final ethanol concentration above 0.1% due to potential cell toxicity.

Published Applications

CANCER RESEARCH

- Suppresses tumor growth in breast cancer cells (Freiss et al.; Osborne et al.).

References

Freiss G et al. (1990) Mechanisms of 4-hydroxytamoxifen anti-growth factor activity in breast cancer cells: alterations of growth factor receptor binding sites and tyrosine kinase activity. *Biochem Biophys Res Commun* 173(3): 919–26.

Osborne CK et al. (1991) Acquired tamoxifen resistance: correlation with reduced breast tumor levels of tamoxifen and isomerization of trans-4-hydroxytamoxifen. *J Natl Cancer Inst* 83(20): 1477–82.

Poulin R et al. (1989) Antiestrogenic properties of keoxifene, trans-4-hydroxytamoxifen, and ICI 164384, a new steroidal antiestrogen, in ZR-75-1 human breast cancer cells. *Breast Cancer Res Treat* 14(1): 65–76.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.