

Small Molecules

MLN4924

Inhibits NEDD8-activating enzyme

Catalog # 74182
74184

1 mg
10 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM


FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

MLN4924 inhibits cullin neddylation by directly targeting the NEDD8-activating enzyme that is essential for the activity of cullin-RING ligase (CRL) (Nawrocki et al.). The components of CRL are frequently up-regulated in human cancers (Zhao & Sun), and MLN4924 was reported to suppress the growth of human tumor xenografts in mice (Soucy et al.).

Molecular Name:	MLN4924
Alternative Names:	Pevonedistat
CAS Number:	905579-51-3
Chemical Formula:	C ₂₁ H ₂₅ N ₅ O ₄ S
Molecular Weight:	443.5 g/mol
Purity:	≥ 98%
Chemical Name:	sulfamic acid, [(1S,2S,4R)-4-[4-[[[(1S)-2,3-dihydro-1H-inden-1-yl]amino]-7H-pyrrolo[2,3-d]pyrimidin-7-yl]-2-hydroxycyclopentyl]methyl ester

Structure:



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 45 mM · Absolute ethanol ≤ 25 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 225 μL of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

GENOME EDITING

· In combination with other small molecules (“CRISPY” mix), MLN4924 has been reported to increase precise CRISPR-Cas9 genome editing in human pluripotent stem cells (Riesenberg & Maricic).

CANCER RESEARCH

· By blocking NAE and proteasome, MLN4924 together with bortezomib inhibits AKT and mTOR and induces apoptosis in multiple myeloma cell lines (Gu et al.).

References

Gu Y et al. (2014) MLN4924, an NAE inhibitor, suppresses AKT and mTOR signaling via upregulation of REDD1 in human myeloma cells. *Blood* 123(21): 3269–76.

Nawrocki ST et al. (2012) MLN4924: a novel first-in-class inhibitor of NEDD8-activating enzyme for cancer therapy. *Expert Opin Investig Drugs* 21(10): 1563–73.

Riesenberg S & Maricic T. (2018) Targeting repair pathways with small molecules increases precise genome editing in pluripotent stem cells. *Nat Commun* 9(1): 2164.

Soucy TA et al. (2009) An inhibitor of NEDD8-activating enzyme as a new approach to treat cancer. *Nature* 458(7239): 732–6.

Zhao & Sun. (2013) Cullin-RING ligases as attractive anti-cancer targets. *Curr Pharm Des* 19(18): 3215–25.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.