

Small Molecules

SMER28

Enhances mammalian autophagy

Catalog # 74202
74204

5 mg
10 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

SMER28, a bromo-substituted quinazoline, was first identified as a small-molecule enhancer of rapamycin (SMER; Sarkar et al.). SMER28 induces autophagy by increasing autophagosome synthesis and enhancing the clearance of model autophagy substrates (Renna et al.).

Molecular Name:	SMER28
Alternative Names:	Not applicable
CAS Number:	307538-42-7
Chemical Formula:	C ₁₁ H ₁₀ BrN ₃
Molecular Weight:	264.1 g/mol
Purity:	≥ 98%
Chemical Name:	6-bromo-N-prop-2-enylquinazolin-4-amine
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 110 mM · Absolute ethanol ≤ 110 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 379 µL of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

REPROGRAMMING

- As a component of a neural reprogramming cocktail, SMER28 reprograms mouse fibroblasts into induced neural stem cell-like cells (Zhang et al.).

DISEASE MODELING

- Enhances clearance of autophagy substrates including mutant huntingtin in Huntington's disease and A53T α -synuclein in familial Parkinson's disease (Sarkar et al.).
- Through autophagy factor ATG5, stimulates erythropoiesis and up-regulates expression of globin genes in both in vitro and in vivo models of Diamond-Blackfan anemia (Doulatov et al.).
- Decreases the level of indicators of Alzheimer's disease (amyloid- β peptide and the amyloid precursor protein [APP]-derived fragment) via the autophagy-related protein 5-dependent pathway (Tian et al.).

References

Doulatov S et al. (2017) Drug discovery for Diamond-Blackfan anemia using reprogrammed hematopoietic progenitors. *Sci Transl Med* 9(376): eaah5645.

Renna M et al. (2010) Chemical inducers of autophagy that enhance the clearance of mutant proteins in neurodegenerative diseases. *J Biol Chem* 285(15): 11061–7.

Sarkar S et al. (2007) Small molecules enhance autophagy and reduce toxicity in Huntington's disease models. *Nat Chem Biol* 3(6): 331–8.

Tian Y et al. (2011) A small-molecule enhancer of autophagy decreases levels of A β and APP-CTF via Atg5-dependent autophagy pathway. *FASEB J* 25(6): 1934–42.

Zhang M et al. (2016) Pharmacological reprogramming of fibroblasts into neural stem cells by signaling-directed transcriptional activation. *Cell Stem Cell* 18(5): 653–67.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.