



EasySep™ Non-Human Primate CD4+ T Cell Isolation Kit

Negative Selection
Catalog #19582

For processing 1 x 10⁹ cells



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713
INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM
FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Document #DX23207 | Version 1_0_0

Description

Isolate untouched and highly purified CD4+ T cells from fresh or previously frozen non-human primate peripheral blood mononuclear cells (PBMCs) in as little as 25 minutes by immunomagnetic negative selection.

- Fast, easy-to-use and column-free
- Up to 85% purity with high recovery
- Untouched, viable cells

This kit targets non-CD4+ T cells for removal with antibodies recognizing specific cell surface markers. Unwanted cells are labeled with antibodies and magnetic particles, and separated without columns using an EasySep™ magnet. Desired cells are simply poured off into a new tube. Isolated cells are immediately available for downstream applications such as flow cytometry, cell culture, or DNA/RNA extraction.

Component Descriptions

| COMPONENT NAME | COMPONENT # | QUANTITY | STORAGE | SHELF LIFE | FORMAT |
|---|-------------|----------|-------------------------------------|--|---|
| EasySep™ Non-Human Primate CD4+ T Cell Isolation Cocktail | 19582C | 1 x 1 mL | Store at 2 - 8°C. Do not freeze. | Stable until expiry date (EXP) on label. | A combination of monoclonal antibodies in PBS and 0.09% sodium azide. |
| EasySep™ D Magnetic Particles | 19250 | 2 x 1 mL | Store at 2 - 8°C. Do not freeze. | Stable until expiry date (EXP) on label. | A suspension of magnetic particles in TBS. |

PBS - phosphate-buffered saline; TBS - Tris-buffered saline

Components may be shipped at room temperature (15 - 25°C) but should be stored as indicated above.

Sample Preparation

This kit has been verified for use with rhesus and cynomolgus macaques.

PERIPHERAL BLOOD

For peripheral whole blood from rhesus macaques, prepare a PBMC suspension by centrifugation over a density gradient medium (e.g. Lymphoprep™, Catalog #07801). For peripheral whole blood from cynomolgus macaques, dilute the density gradient medium to 90% using D-PBS (Without Ca⁺⁺ and Mg⁺⁺; Catalog #37350).

NOTE: For higher recovery, 14 mL polystyrene round-bottom tubes (e.g. Catalog #38008) are recommended for density gradient centrifugation, particularly for smaller volumes of peripheral blood.

For samples > 24 hours old, it may be necessary to lyse the red blood cells (RBCs) using Ammonium Chloride Solution (Catalog #07800) prior to cell isolation.

If using previously frozen PBMCs, incubate the cells with DNase I Solution (Catalog #07900) at a concentration of 100 µg/mL at room temperature (15 - 25°C) for at least 15 minutes prior to labeling and separation. Filter aggregated suspensions through a 37 µm cell strainer (e.g. Catalog #27250) for optimal results.

After preparation, resuspend cells at 5 x 10⁷ cells/mL in recommended medium.



Recommended Medium

EasySep™ Buffer (Catalog #20144), RoboSep™ Buffer (Catalog #20104), or PBS containing 2% FBS and 1 mM EDTA. Medium should be free of Ca⁺⁺ and Mg⁺⁺.

Directions for Use – Manual EasySep™ Protocols

See page 1 for Sample Preparation and Recommended Medium. Refer to Tables 1 and 2 for detailed instructions regarding the EasySep™ procedure for each magnet.

Table 1. EasySep™ Non-Human Primate CD4+ T Cell Isolation Kit Protocol

| | | EASYSEP™ MAGNETS | |
|------|--|--|--|
| STEP | INSTRUCTIONS |  EasySep™ (Catalog #18000) |  “The Big Easy” (Catalog #18001) |
| 1 | Prepare sample at the indicated cell concentration within the volume range. | 5 x 10 ⁷ cells/mL 0.5 - 2 mL | 5 x 10 ⁷ cells/mL 0.5 - 6 mL |
| | Add sample to required tube. | 5 mL (12 x 75 mm) polystyrene round-bottom tube (e.g. Catalog #38007) | 14 mL (17 x 95 mm) polystyrene round-bottom tube (e.g. Catalog #38008) |
| 2 | Add Isolation Cocktail to sample. | 50 µL/mL of sample | 50 µL/mL of sample |
| | Mix and incubate. | RT for 10 minutes | RT for 10 minutes |
| 3 | Vortex EasySep™ D Magnetic Particles. NOTE: Particles should appear evenly dispersed. | 30 seconds | 30 seconds |
| 4 | Add Magnetic Particles to sample. | 100 µL/mL of sample | 100 µL/mL of sample |
| | Mix and incubate. | RT for 5 minutes | RT for 5 minutes |
| 5 | Add recommended medium to top up the sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times. | Top up to 2.5 mL | <ul style="list-style-type: none"> • Top up to 5 mL for samples < 2 mL • Top up to 10 mL for samples ≥ 2 mL |
| | Place the tube (without lid) into the magnet and incubate. NOTE: Purity may be improved by increasing incubation time to 10 minutes; however, this may reduce recovery. | RT for 5 minutes | RT for 5 minutes |
| 6 | Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring the cell suspension into a new tube. | Use a new 5 mL tube | Use a new 14 mL tube |
| 7 | Place the new tube from step 6 (without lid) into the magnet and incubate for a second separation. | RT for 5 minutes | RT for 5 minutes |
| 8 | Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring the cell suspension into a new tube. | Use a new 5 mL tube Isolated cells are ready for use | Use a new 14 mL tube Isolated cells are ready for use |

RT - room temperature (15 - 25°C)

* Leave the magnet and tube inverted for 2 - 3 seconds, then return upright. Do not shake or blot off any drops that may remain hanging from the mouth of the tube.

Directions for Use – Fully Automated RoboSep™ Protocol

See page 1 for Sample Preparation and Recommended Medium. Refer to Table 3 for detailed instructions regarding the RoboSep™ procedure.

Table 3. RoboSep™ Non-Human Primate CD4+ T Cell Isolation Kit Protocol

| STEP | INSTRUCTIONS | RoboSep™ (Catalog #21000) |
|------|--|---|
| 1 | Prepare sample at the indicated cell concentration within the volume range. | 5 x 10 ⁷ cells/mL 0.5 - 6 mL |
| | Add sample to required tube. | 14 mL (17 x 95 mm) polystyrene round-bottom tube (e.g. Catalog #38008) |
| 2 | Select protocol. | Non-Human Primate CD4+ T Cell Isolation 19582 |
| 3 | Vortex EasySep™ D Magnetic Particles. NOTE: Particles should appear evenly dispersed. | 30 seconds |
| 4 | Load the carousel. | Follow on-screen prompts |
| | Start the protocol. | Press the green "Run" button |
| 5 | Unload the carousel when the run is complete. Remove the tube containing the isolated cells. | Isolated cells are ready for use |

Notes and Tips

ASSESSING PURITY


NOTE: Due to the presence of residual RBCs, use of a CD45 antibody is strongly recommended.

NOTE: Use of a cell viability dye is strongly recommended.

For purity assessment of non-human primate CD4+ T cells (CD3+CD4+) by flow cytometry, use the following fluorochrome-conjugated antibodies:

- Anti-human CD3 antibody, clone SP34.2,
- Anti-Human CD4 Antibody, Clone OKT4 (Catalog #60016),
- Anti-Human CD8 Antibody, Clone SK1 (Catalog #60022), and
- Anti-human CD45 antibody, clone D058-1283

Data



In the above example, starting with rhesus macaque PBMCs, the CD4+ T cell content (within the CD45+CD3+ gate) of start and final isolated fractions are 49.5% and 92.3%, respectively.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, RoboSep, and SepMate are trademarks of STEMCELL Technologies Canada Inc. Lymphoprep is a trademark of Alere Technologies. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.