

# Axygen® PCR Plasticware Compatibility Chart



PCR/qPCR THERMOCYCLER		AXYGEN® PCR PLASTICWARE									
		MICROPLATE								TUBE	
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL
		CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094
<b>Agilent®</b>	SureCycler® 8800	•	•	•	•	•	•	•	•	•	•
<b>Analytik Jena/ Biometra</b>	FlexCycler	•	•	•	•	•	•	•	•	•	•
	T1 Thermal Cycler	•	•	•	•	•	•	•	•	•	•
	TProfessional/Basic	•	•	•	•	•	•	•	•	•	•
	TProfessional/ Standard	•	•	•	•	•	•	•	•	•	•
	TProfessional/TRIO	•	•	•	•	•	•	•	•	•	•
	TRobot	•	•	•	•	•	•	•	•	•	•
	Uno II	•	•	•	•	•	•	•	•	•	•
	TGradient	•	•	•	•	•	•			•	•
	TPersonal	•	•	•	•	•	•			•	•
	Uno	•	•	•	•	•	•			•	•
	qTOWER	•	•	•	•		•			•	•
	SpeedCycler <sup>2</sup>	•	•	•	•		•			•	•
Optical	•	•	•	•		•			•	•	

		AXYGEN® PCR PLASTICWARE									
		MICROPLATE								TUBE	
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL
PCR/qPCR THERMOCYCLER		CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094
Applied Biosystems®	9700	•		•	•		•	•	•	•	•
	Fast Veriti®/Fast	•		•	•	•*	•	•	•	•	•
	7900HT/Fast			•	•	•*	•		•	•	•
	ViiA™ 7/Fast			•	•	•*	•		•	•	•
	2720	•		•	•		•			•	•
	9600	•		•	•		•			•	•
	7500/Fast			•	•	•*	•			•	•
	7000			•	•		•			•	•
	7300			•	•		•			•	•
	7700			•	•		•			•	•
	StepOnePlus™					•			•	•	•
	ProFlex™ PCR System			•	•		•			•	•
	QuantStudio™ 12K			•	•		•			•	•

\*Compatible with ABI® Fast PCR only.

PCR/qPCR THERMOCYCLER		AXYGEN® PCR PLASTICWARE											
		MICROPLATE								TUBE			
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL		
CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094				
Applied Biosystems® (continued)	9800 Fast												
	StepOne™												
Axygen®	MaxyGene™ II		•	•	•	•	•				•	•	
BioRad®/MJ Research®	C1000 Touch™	•	•	•	•		•	•	•	•	•	•	•
	DNA Engine Dyad®/ Dyad Disciple™	•	•	•	•		•	•	•	•	•	•	•
	Engine Tetrad® 2	•	•	•	•		•	•	•	•	•	•	•
	PTC-100®	•	•	•	•		•	•	•	•	•	•	•
	PTC-200	•	•	•	•		•	•	•	•	•	•	•
	PTC-225 Tetrad	•	•	•	•		•	•	•	•	•	•	•
	S1000™	•	•	•	•		•	•	•	•	•	•	•
	T100™	•	•	•	•		•	•	•	•	•	•	•
	MyCycler™	•	•	•	•		•				•	•	
	CFX Connect™	•	•	•	•	•	•				•	•	
CFX96™ Touch	•	•	•	•	•	•				•	•		

		AXYGEN® PCR PLASTICWARE									
		MICROPLATE								TUBE	
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL
PCR/qPCR THERMOCYCLER		CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094
<b>BioRad®/MJ Research®</b> (continued)	Chromo4™	•	•	•	•		•			•	•
	DNA Engine Opticon® 2	•	•	•	•		•			•	•
	iCycler	•	•	•	•		•			•	•
	iQ™5	•	•	•	•		•			•	•
	MiniOpticon™	•	•	•	•		•			•	•
	MyiQ™	•	•	•	•		•			•	•
	CFX384 Touch™							•	•		
	Opticon™		•			•					
	Opticon™ 2		•			•					
<b>Corbett™ Life Science</b>	Palm-Cycler™ 96	•	•	•	•		•			•	•
	Palm-Cycler™ 384							•	•		

PCR/qPCR THERMOCYCLER		AXYGEN® PCR PLASTICWARE									
		MICROPLATE								TUBE	
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL
		CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094
Eppendorf®	Mastercycler® pro/S	•	•	•	•		•	•	•	•	•
	Mastercycler®	•	•	•	•		•			•	•
	Mastercycler® ep	•	•	•	•		•			•	•
	Mastercycler® nexus	•	•	•	•		•			•	•
	Mastercycler® nexus eco	•	•	•	•		•			•	•
	Mastercycler® nexus gradient	•	•	•	•		•			•	•
	Mastercycler® nexus gradient eco	•	•	•	•		•			•	•
	Mastercycler® nexus GSX1	•	•	•	•		•			•	•
	Mastercycler® nexus GSX1e	•	•	•	•		•			•	•
	Mastercycler® nexus SX1	•	•	•	•		•			•	•
	Mastercycler® nexus SX1e	•	•	•	•		•			•	•

PCR/qPCR THERMOCYCLER		AXYGEN® PCR PLASTICWARE									
		MICROPLATE								TUBE	
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL
CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094		
<b>Eppendorf®</b> (continued)	Mastercycler® ep realplex	•	•	•	•		•		•	•	
	Mastercycler® M384							•	•		
	Mastercycler® pro 384							•	•		
	Mastercycler® nexus flat								•		
	Mastercycler® nexus flat eco								•	•	
<b>Ericom</b>	Deltacycler I	•		•	•	•	•		•	•	
	Single Block	•		•	•	•	•		•	•	
	Twin Block	•		•	•	•	•		•	•	
<b>MWG™</b>	Primus 96	•	•	•	•	•	•		•	•	
	Primus 384							•	•		
<b>PeqLab™</b>	peqSTAR 96		•	•	•		•		•	•	
	peqSTAR 96Q		•	•	•		•		•	•	
	peqSTAR 384							•	•		

PCR/qPCR THERMOCYCLER		AXYGEN® PCR PLASTICWARE									
		MICROPLATE								TUBE	
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL
CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094		
Qiagen®	FlexiGene®	•	•	•	•		•	•	•	•	•
	Rotor-Gene® Q								•		
Roche	Lightcycler® 480										
SensoQuest™	Thermocycler 384		•	•	•		•	•	•	•	•
	Thermocycler 96		•	•	•		•	•	•	•	•
Stratagene™	Robocycler®	•	•	•	•		•		•	•	
	MX3000®			•	•		•		•	•	
	MX3000P®			•	•		•		•	•	
	MX3005P®			•	•		•		•	•	
	MX4000®			•	•		•		•	•	
TaKaRa™	TP3000	•	•	•	•	•	•		•	•	
	TP240		•	•	•		•		•	•	

PCR/qPCR THERMOCYCLER		AXYGEN® PCR PLASTICWARE									
		MICROPLATE								TUBE	
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL
		CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094
Technē™	Touchgene™ Gradient (TC512)	•	•	•	•	•	•	•	•	•	•
	Touchgene™ X	•	•	•	•	•	•	•	•	•	•
	Genius		•	•	•		•	•	•	•	•
	Genius Quad		•	•	•		•	•	•	•	•
	Genius (TC412)	•	•	•	•	•	•			•	•
	Prime Elite		•	•	•		•	•	•	•	•
	Prime Elite Satellite		•	•	•		•	•	•	•	•
	PrimeG		•	•	•		•	•	•	•	•
	Cyclogene™		•	•	•		•			•	•
	Flexigene™		•	•	•		•			•	•
	Prime		•	•	•		•	•	•	•	•
	PrimeQ		•	•	•		•			•	•



		AXYGEN® PCR PLASTICWARE									
		MICROPLATE								TUBE	
		96-well, No Skirt	96-well, Full Skirt	96-well, Half Skirt, Double Notch	96-well, Flat Top, No Skirt	96-well, Flat Top, Low Profile, No Skirt	96-well, Half Skirt, Single Notch	384-well, Full Skirt, Double Notch	384-well, Full Skirt, ABI Compatible	PCR Tube with Flat Cap, 0.2 mL	PCR Strip Tubes (8/Strip) 0.2 mL
PCR/qPCR THERMOCYCLER		CATALOG # 38098	CATALOG # 38099	CATALOG # 38100	CATALOG # 38101	CATALOG # 38102	CATALOG # 38103	CATALOG # 38105	CATALOG # 38106	CATALOG # 38091	CATALOG # 38094
<b>Techne®</b> (continued)	Quantica®		•	•	•			•		•	•
	3Prime									•	•
	3PrimeG									•	•
	3PrimeX									•	•
<b>Thermo Hyaid</b>	Multiblock System	•	•	•	•	•	•	•	•	•	•
	Omnigene	•	•	•	•	•	•	•	•	•	•
	PCR Express	•	•	•	•	•	•	•	•	•	•
	Px2	•	•	•	•	•	•	•	•	•	•
	PxE	•	•	•	•	•	•	•	•	•	•
	Touchdown	•	•	•	•	•	•	•	•	•	•
	Omn-E	•	•	•	•	•	•			•	•
	PCR Sprint									•	•

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.